

ThermoFisher
SCIENTIFIC

Laboratory Solutions India

**Best Performance - 2020
Consumables – Platinum Category**

Shan Biotech & Diagnostics

4th February 2021, Mumbai





**ART tips are now available
for ALL your applications:**


- BARRIER
- NON-FILTERED
- SPECIALTY



Finnpipette F2



Excellence in Form and Function



Safe and Comfortable

Novus excels in repetitive applications – such as plate filling, due to its multidispensing functionality. Electronic pipetting operation and adjustable speed make it ideal when working with small volumes and viscous liquids.

The combination of its lightweight design and soft-touch tip ejection makes Novus the safe and comfortable choice for heavy-duty pipetting applications. Its unique index-finger operation allows your thumb to relax during pipetting. Novus also ensures your data integrity by storing programs and monitoring the number of pipetting steps between calibration intervals to provide quality, consistent results.



Lightweight handle design combined with low plunger tip attachment and ejection forces help reduce the risk of RSI.

Adjustable comfort

The 120° adjustable finger rest – for both right- and left-handed operator comfort – allows your hand to relax during pipetting cycles.

Comfortable Tip Ejection

Soft touch tip ejection combined with the softer collar of Finntip Flex tips reduces ejection forces.

Super blow-out

150% increase in air boost ensures efficient delivery of micro-volumes and prevents capillary action in 50 µL models and below.

NEW silent, effortless volume adjustment

Enjoy the new silent and light volume adjustment with over 50% lower forces compared to earlier models. Improved grip on the volume adjustment knob enables effortless and precise adjustment.

Secure volume lock

Securely set and lock volumes to prevent undesirable volume drift during operation.